



The HDice Target at JLab



Collaborators

Jefferson Lab (HDice)

G. Dezern, C. Hanretty, T. Kageya, M. Lowry, A. M. Sandorfi, X. Wei

Jefferson Lab (UITF)

H. Areti, J. Grames, J. Gubeli, M. Poelker, W. Akers

<u>Universita di Roma Tor Vergata and INFN-Sezione di Roma2</u> (Gas analysis)

A. D'Angelo

Universita di Ferrara and INFN di Ferrara (MgB2 magnet)

M. Contalbrigo, P. Lenisa, M. Statera




Overview

- Uses of HDice target
- Production of a target cell for experiment
- Polarization process/mechanisms
- Lessons learned: eHD during g14
- eHD test in the Upgraded Injector Test Facility



HDice


- A new type frozen spin target.
- Target material consists solely of polarizable protons and neutrons ⇒ no dilution factor coming from target material.
- Target material possesses a T1 on the order of years (no repolarization needed).
- Is a very complicated system requiring many steps in the production of a polarized target.






What is it good for?

- Has been used (with photons) in Hall B as part of the N* program
 ⇒ g14 (Nov 2011-May 2012)
- N* program seeks/sought to explore the excited states of the nucleon to better understand the internal dynamics of protons and neutrons





$$I = I_0 \left\{ \left(1 + \overrightarrow{\Lambda}_i \cdot \overrightarrow{\mathbf{P}} \right) + \delta_{\odot} \left(\mathbf{I}^{\odot} + \overrightarrow{\Lambda}_i \cdot \overrightarrow{\mathbf{P}^{\odot}} \right) + \delta_I \left[\sin(2\beta) \left(\mathbf{I}^{\mathbf{s}} + \overrightarrow{\Lambda}_i \cdot \overrightarrow{\mathbf{P}^{\mathbf{s}}} \right) + \cos(2\beta) \left(\mathbf{I}^{\mathbf{c}} + \overrightarrow{\Lambda}_i \cdot \overrightarrow{\mathbf{P}^{\mathbf{c}}} \right) \right] \right\}$$




What is it good for?

- Next up: Transversely polarized frozen spin target for use with electrons.
- Three (A-rated) proposals from PAC 41 rated as having a high impact for Hall B:
 - SIDIS, C12-11-111, Marco Contalbrigo,... ⇒ [A;C1]
 - Dihadron production, PR-12-009, Harut Avakian,... ⇒ [A;C1]
 - DVCS, PR12-12-101, Latifa Elouadrhiri,...
 ⇒ [A;C1]
- C1

 ⇒ requires successful demonstration to Lab management of viable performance in a subsequent eHD test run



 Provides a way to study proton tomography through the measurement of Generalized Parton Distributions (GPDs)






 Provides a way to study proton tomography through the measurement of Generalized Parton Distributions (GPDs):

$$E(x,\xi,t)$$
, $\tilde{E}(x,\xi,t)$, $H(x,\xi,t)$, $\tilde{H}(x,\xi,t)$


- GPDs themselves provide access to the underlying framework; observables are calculated via integration of GPDs
- Tomographic distributions:



Transverse target polarization asymmetries are **required** for access to **E**




(Why not a conventional target?)





(Why not a conventional target?)




No large magnets to limit acceptance!




Life Cycle of HDice Target

Gas Handling




PicHDice polarized target production

- Made in the Production Dewar (PD)
 - gas to solid
 - initial NMR (TE calib)
- Transferred to Oxford Dilution Fridge (DF)
 - Polarized in a 15 T B field at 10 mK for at least 3 months
- Transferred to PD for NMR measurements
- Transferred to Storage Dewar (SD) for storage (or back to the DF)
- Transferred to the In-Beam Cryostat (IBC) for data taking
- All transfers facilitated by Transfer Cryostat (TC)






Target Transfers









Target transfer between PD and DF for polarization





Target transfer between PD and DF for polarization








Target transfer between PD and DF for polarization



TC mounted on top of the DF











How does one make polarized, frozen-spin HD?






What polarization mechanisms can we use?

- For all molecules, we have 2 electrons paired in 1s orbital ⇒ no DNP
- For all molecules, we (initially) have $S \neq 0 \Rightarrow$ couples to B field






External magnetic field rapidly aligns the $Ortho-H_2$ and $Para-D_2$ which spin-exchange with the H and D (respectively) in the HD crystal.








This process continues until there are only inactive, L=0 states remaining.

The nucleons are now in a frozen-spin state: a state in which the polarized nucleons are so cold, there are (essentially) no available depolarization mechanisms.









 H_{2}


HD

 D_{2}




(in words)

- Solid HD target is loaded into Dilution Fridge.
- Spins are aligned using a 15 T magnet at temperatures of ~10 mK.
- The (L=1) H₂ and D₂ molecules are polarized (are considered impurities).
- Polarized (L=1) H₂ and D₂ nuclei spin exchange with HD, polarizing it.
- After ~3 months, the L=1 H₂ and D₂ molecules decay down to L=0 and are no longer polarizable.
- TE for these conditions lead to P(H) = 90% and P(D) = 30%






Polarizing the Target






Polarizing the Target

- After 3 months, target transferred from DF to PD for NMR
- Polarization can be determined and compared to before HD was aged in DF




Charles Hanretty (JLab)






Target transfer to SD for storage/transport






Target transfer to SD for storage/transport





Target transfer to SD for storage/transport






Target transfer to SD for storage/transport



Target Transport

 SD (with target) may be transported while magnet is on (persistent mode) and LHe space pumped to 1.7 K






Transporting SD to Hall B (g14)

Crane lift of SD with target inside (Hall B)



In-Beam Cryostat (IBC)


 For in-beam operations (data-taking), the In-Beam Cryostat (IBC) is used.






In-Beam Cryostat (IBC)


- IBC is a Dilution Refrigerator capable of operating both vertically (for docking with TC) and horizontally (for data-taking).
- T = 50 mK
- B = 1 T (solenoid); 0.075 T (saddle)









In-Beam Cryostat (IBC)



T = 50 mK

B = 1 T (solenoid); 0.075 T (saddle)







Ready for running but...

- ...Is there anything we can do to improve the target polarization?
- P(H) = 90% and P(D) = 30% are <u>ideal</u> values.
- In reality, $P(H) = 60\pm5\%$ and $P(D) = 15\pm7\%$ is essentially the best we can hope for using aging alone.
- Can we give some of the H polarization to D?




Polarization Transfer

- Forbidden Adiabatic Fast Passage (FAFP):
 - Uses one large RF field sweep to transfer P(H) to P(D)
 - Pros: potential 100% efficiency
 - Cons: only one chance (one large RF field sweep), everything must work perfectly, high RF power heats fridge, (RF) field jitter decreases efficiency
- Saturated Fast Passage-Multi FAFP:
 - Uses many small RF sweeps to transfer P(H) to P(D)
 - Pros: guaranteed efficiency is 50%, lower RF power → less heating,
 (RF) field jitter helps
 - Cons: max efficiency is 50%



Polarization Transfer





Analyses using g14 data

Dao Ho (CMU):

$$\gamma_c D \rightarrow p \pi^-(p)$$

$$\gamma_c D \to K^0 \Lambda (p)$$

$$\gamma_c D \rightarrow D \pi^+ \pi^-$$

$$\gamma_c D \to K^0 \Sigma^0$$

Irene Zonta (Roma-II):

$$\gamma_c n(p) \rightarrow \rho n(p) \rightarrow n \pi^+ \pi^-(p)$$

Haiyun Lu (CMU, U. Iowa):

$$\gamma_L n(p) \rightarrow p \pi^-(p)$$

Peng Peng (UVA):

$$\gamma_c p(n) \rightarrow p \pi^+ \pi^-(n)$$

$$\gamma_c n(p) \rightarrow n \pi^+ \pi^-(p)$$

$$\gamma_c n(p) \rightarrow p \pi^-(p)$$

Tsuneo Kageya (JLab):

$$\gamma_c n(p) \rightarrow p \pi^-(p)$$

Jamie Fleming (U. Edinburgh):


$$\gamma_{\tau} p \rightarrow p \pi^{+} \pi^{-}$$

$$\gamma_L n(p) \rightarrow n \pi^+ \pi^-(p)$$

$$\gamma_{c}n(p) \rightarrow K^{+} \Sigma^{-}(p)$$



Analyses using g14 data



Peng Peng (UVA):

$$\gamma_c p(n) \rightarrow p \pi^+ \pi^-(n)$$

$$\gamma_c n(p) \rightarrow n \pi^+ \pi^-(p)$$

$$\gamma_c n(p) \rightarrow p \pi^-(p)$$

Tsuneo Kageya (JLab):

To be discussed at the Hadron Spectroscopy Working Group Meeting tomorrow


$$\gamma_L n(p) \rightarrow p \pi^-(p)$$

$$\gamma_L p \rightarrow p \pi^+ \pi^ \gamma_L n(p) \rightarrow n \pi^+ \pi^-(p)$$
 $\gamma_C n(p) \rightarrow K^+ \Sigma^-(p)$



eHD

 During g14, eHD tests were conducted to study how the target (and IBC) fared with electrons.






eHD during g14

- Three depolarization mechanisms were found.
- Local beam heating:
 - The T1 of the target material is a function of B/T.
 - Temperature spikes reduce the $T1 \rightarrow$ no more frozen spin
- Unpairing of 1s electrons:
 - Beam heating unpairs the electrons in 1s shell
 - Residual electron is unpolarized (function of temp) → flips with a frequency that has Fourier components at D and H Larmor frequencies
 - Depolarizes nearby HD; depolarization spreads
- Ionization of HD:
 - HD molecule loses an electron (becomes ionized), is highly reactive
 - Chain of chemical reactions begin which culminate in a "recombination flash" → produces significant heat in target → loss of polarization



eHD during g14





eHD during g14

- Three depolarization mechanisms were found.
- Local beam heating:
 - The T
 - Temp
- Unpairin
 - Beam
 - Residfrequent

Solutions:

- Make target cell shorter and fatter
- Make Aluminum wires shorter and fatter to remove heat faster
- Faster Raster

yuencies

a

- Depolarizes nearby HD; depolarization spreads
- Ionization of HD:
 - HD molecule loses an electron (becomes ionized), is highly reactive
 - Chain of chemical reactions begin which culminate in a "recombination flash" → produces significant heat in target → loss of polarization



eHD: Target cell redesign


g14 Cell:





eHD: Target cell redesign


- eHD target cell will be half the length and almost twice the diameter of the g14 cell
- Aluminum wires will be shorter and thicker
- Effect: Lower equilibrium temperature of HD






eHD: Target cell redesign

- eHD target cell will be half the length and almost twice the diameter of the g14 cell
- Aluminum wires will be shorter and thicker
- Effect: Lower equilibrium temperature of HD





eHD: Faster Raster



Charles Hanretty (JLab)


The HDice Target at JLab

17 June 2015



eHD: Faster Raster

 New target cell design leads to equilibrium temp of 145 mK but local temps can still reach over 1000 mK




Charles Hanretty (JLab) The HDice Target at JLab 17 June 2015



eHD: Faster Raster


- New target cell design leads to equilibrium temp of 145 mK but local temps can still reach over 1000 mK
- By increasing the speed of the Raster (and altering the pattern),
 the temperature of the target material is greatly reduced



Charles Hanretty (JLab) The HDice Target at JLab 17 June 2015



eHD Tests in UITF




Charles Hanretty (JLab)

The HDice Target at JLab

17 June 2015



eHD Tests in UITF




Charles Hanretty (JLab)

The HDice Target at JLab

17 June 201



eHD Tests in UITF





Summary

- HDice is a new type frozen spin target designed for fixed target experiments at JLab.
- Is a complicated target system utilizing many steps and pieces of equipment in the production of a single target cell.
- Has been used with photons as part of the N* program in Hall B (g14).
- Has potential as a transversely polarized frozen spin target for use with electrons → no large magnets.
- eHD tests to be carried out in the UITF in Fall 2016




END



Cryostats


Helium Evaporation Fridge





Cryostats

(Horizontal) Dilution Fridge





Great, so that means that P(H)=90% and P(D)=30%, right?

If it were only that simple

P(H)=90% and P(D)=30% are **TE** values at 15 T and 10 mK

There are two effects that must strike a balance in this process:

As the (L=1) H2 and D2 states decay, the T1 of the target material grows


⇒ takes longer to reach the TE values (deg of pol increases slower and slower) L=1→0 decays release heat which must be removed from the HD in order to reach the TE values

⇒ HD can't immediately reach low temperatures (10 mK) (max pol is decreased)

The result is that $P(H)=60\pm5\%$ and $P(D)=15\pm7\%$




Target transfer & transport




Target transfer to SD for storage/transport



TC mounted on top of the DF







Charles Hanretty (JLab) The HDice Target at JLab 17 June 2015



Target stacking in DF








UITF Beamline


(A compressed view)






UITF Beamline

(A compressed view)





The following slides are pictures/cartoons made by [CMH]















 $H_{_{2}}$




HD




 D_2












Cryostats

Helium Evaporation Fridge





Cryostats

(Horizontal) Dilution Fridge

